

Catheter Valves

This product should not be used without assessment of bladder function.

Contra-indications:

- Reduced bladder capacity
- No bladder sensation
- Cognitive impairment
- Insufficient manual dexterity to operate the catheter valve

Valve should be in place for 5-7 days.

Prescribing Advice: 1 box of 5 should be sufficient for 1 month supply.

Unomedical Ltd

Careline Lever catheter valve

90 821 1S-20

Penile sheaths

These must be measured and fitted by appropriately trained staff and measuring should only be with manufacturer's measuring guide.

Prescribing Advice: 1 box of 30 should be sufficient for 1 month supply.

Manfred Sauer UK Ltd

P-SURE Latex free self-adhesive

Paediatric	18mm	97.18	Small -	20mm	97.20
Small	22mm	97.22	Small+	24mm	97.24
Medium -	26mm	97.26	Medium	28mm	97.28
Medium +	30mm	97.30	Large -	32mm	97.32
Large	35mm	97.35	Large +	37mm	97.37
Extra Large	40mm	97.40			

Rochester Medical Ltd

Freedom Transfix Silicone self adhesive

Standard	25mm	TF12530	29mm	TF12930
	32mm	TF13230	36mm	TF13630
	41mm	TF14130		
Short	25mm	TF22530	29mm	TF22930
	32mm	TF23230	36mm	TF23630
	41mm	TF24130		
Extra wide adhesive band	25mm	TF32530	29mm	TF32930
	32mm	TF33230	36mm	TF33630
	41mm	TF34130		

Urinary Catheter Formulary



Indwelling Catheters

Pre-filled catheters are not recommended for use and any patients currently prescribed these should be switched to one of the foley catheters in the formulary.

Catheter trays should no longer be prescribed or used as they are not cost effective.

Short/medium term use (1-3 weeks) and no latex allergy

Rusch PTFE Aquaflate (Teleflex Medical) sizes 12-24 10 ml balloon
Coated latex with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation.

Male	DP3101 12-24
Female	DP2101 12-24

Long term use (3-12weeks)

Prescribing advice: up to 3 catheters can be prescribed at one time; they should be prescribed as required, about every 3- 8 months, if required more frequently an assessment may be required.

Medasil Surgical Ltd All Silicone

All silicone catheter, sizes 12-26, 10 ml balloon, not licensed for suprapubic use.

Male	84M	12-26
Female	86F	12-26

Suprapubic catheters

L.IN.C Medical Systems Ltd

2 way all silicone suprapubic catheter with integral balloon and central opening, sizes 16-24 with 10ml balloon.

Short	0847 1610-2410
Long	0845 1610-2410

Where silicone balloon is causing problems with removal and there is no latex allergy

Rusch PTFE SympaCath AquaFlate (Teleflex Medical) sizes 16-24 with 10ml balloon with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation. This catheter is not open ended.

Female (23cm length)	DH2101 16-24
Standard (41cm length)	DH3101 16-24

Lidocaine and Chlorhexadine Gel

Please be aware of the contraindications for use and observations required for the use of lidocaine. It needs to be left in place for 4-5 minutes to be effective.

Catheter Leg Drainage Bags

Bags should be used for 5-7 days

Clinisupplies Ltd

ProSys Sterile leg bags contain 1 pair of non-sterile latex gloves and soft elasticised cotton Velcro straps per box 10.

Prescribing advice: 1 box of 10 should be prescribed every 2 months only.

Slide Action Tap

Short tube	350ml	P350S
	500ml	P500S
	750ml	P750S
Long tube	350ml	P350L
	500ml	P500L
	750ml	P750L

Lever Tap

Short tube	500ml	P500S-LT
Long tube	500ml	P500L-LT

Catheter Night Drainage bags

Clinisupplies Ltd

ProSys Sterile 2 litre bag with 90cm inlet tube, non-return valve, sample port connector, slide outlet tap and 1pair non latex gloves, 10 per box.

Prescribing advice: 1 box of 10 bags should only be prescribed every 2 months; bags are reusable for 5-7days, if more required then single use bags may be more cost effective

ProSys	P2000
ProSys with lever tap outlet	P2000-LT

Clinisupplies Ltd

Non-sterile 2 litre night bag with single use drainable tap, single use, discard following emptying, 10 per box

Prescribing advice: 3 boxes of 10 bags should be sufficient for 1 month supply.

ProSys	PSU2
--------	------